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CuroCell® 4

Examination of the pressure relief capability and the shear force

(Berlin Cert, 02/09/09)

CuroCell® 4 gets highest rate, "High Pressure Relief", during testing of the product of accredited institutions (29 % better compared to the reference mattress at weight of 80 kg).

Pressure Relief Index (PRI) (Load 80 kg)

100 % below 30 mmHg, 44 % below 20 mmHg, 15 % under 10 mmHg

Pressure Area Index (PAI) (Load 80 kg)

100 % below 30 mmHg, 91 % below 20 mmHg, 47 % below 10 mmHg

Evaluation of the effectiveness of the system

(Care of Sweden, 30/09/09)




Study aimed to describe the function and performance of alternating pressure mattresses feature. This was conducted with test subjects in a weight range between 57 kg to 135 kg. Below are the pressure relief periods presented, ie percentage of a cycle period where the pressure is below threshold values.

Pressure relief period

75 % below 30 mmHg 66 % below 20 mmHg 45 % below 10 mmHg

Performance example

The tables show a cycle period at different loads.

Case Study

History

Female, born 1943 with Multiple Sclerosis (MS). Point according to the Modified Norton scale is 19, BMI is 22. Pressure ulcer category 4 in the sacrum, R/T poor general condition after a respiratory infection. Incontinent both of urine and faeces uses diapers. She is sitting up during the day in a wheelchair, but rests about two hours each day.

Action

She was ordained an alternating pressure relief mattress, CuroCell 4®. She also had a pressure relieving cushion in her wheelchair. Has home-care assistance for position changes at nights.

07/09/07

4 weeks later

11 weeks later

18 weeks later